



DATA SHEET

TISBE Tissue Processor

COMPACT TISSUE PROCESSOR

MODEL: Vacuum tissue processor

CODE: 539001





DATA SHEET

TISBE Tissue Processor

IMAGES:

- 1) OVERVIEW
- 2) PROCESSING CHAMBER
- 3) SOFTWARE

FEATURES:

Processor for histologic samples

Equipped with 12 tanks, 2.5 l with quick connectors

- ✓ 7 reagents tanks
- ✓ 2 cleaning reagents tanks
- ✓ 3 paraffin tanks
- Processing chamber completely stainless steel made and hermetic seal;
- Vacuum and pressure system during processing cycles;
- Processing capacity 300 cassettes;
- Reagent control and management by software;
- Running ambient temperature from 10° to 35° C adjustable;
- Wax working temperature from 52° to 65° C adjustable;
- Visual alarm of reagent saturation;
- Disposable reagents and paraffin tanks with quick connectors for an easy reagents replacement;
- Access to functions through graphical icons on LCD touch screen;



DATA SHEET

TISBE Tissue Processor

- Equipped with two baskets of 150 cassettes;
- 12 different storable processing programs;
- programming in English;
- easy set up of processing times, pressure and vacuum;
- adjustable agitation on each program step;
- personalization of protocols and reagents names;
- suction and charcoal filtering device;
- alarm of charcoal filter saturation;
- sample rescue in case of black out;
- remote alarm: phone warning controller (optional).

DIMENSIONS:

MEASURES W x D x H 1050 x 600 x 580 mm

WEIGHT:

about 60 Kg

TECHNICAL DATA:

Voltage: 110 / 230 V

Frequency: 50 / 60 Hz

Power: 1000 W

Temperature control precision: +/- 1°C

IEC1010 classification: class 1

Instr. Manual: in English

Diapath s.r.l.

Via Salvoldini, 71 24057 Martinengo (BG) Italy
Tel. (+39) 0363.986.411 - Fax (+39) 0363.948.000
www.diapath.com - info@diapath.com



DATA SHEET

TISBE Tissue Processor

EQUIPPED ACCESSORIES:

- | | |
|---------------------------------|-------------|
| - n. 12 paraffin/reagents tanks | Cod. 539010 |
| - n. 2 baskets of 150 cassettes | Cod. 539011 |
| - n. 2 Charcoal filters | Cod. 539012 |
| - n. 1 Basket carrier | Cod. 539013 |

OPTIONAL:

- | | |
|--|-------------|
| - Onda 20 UPS | Cod. 539014 |
| - Uninterruptible power supplier managed by a software | Cod. 539015 |
| - External printer | Cod. 539016 |
| - External fume duct | Cod. 539017 |
| - Auto Dialer | Cod. 539018 |
| - External keyboard | Cod. 539019 |
| - Thermowax paraffin dispenser | Cod. 539020 |
| - Basket bulk (250 Cassettes) | Cod. 539021 |