



CRYOALFA

CRYOSURGERY IN COAT  
POCKET SIZE



- Simple handling
- High efficiency and effectiveness
- Reproducible treatment effects
- Excellent cost/performance ratio

**Cryoalfa Europe Ltd**  
Bernhard-Voß-Str. 25 A  
01445 Radebeul

**Cryoalfa Europe Ltd**

Bernhard-Voß-Str. 25 A  
01445 Radebeul

Tel.: +49 351 795 1866


Fax: +49 351 795 1803

Mail: [medica@cryoswiss.com](mailto:medica@cryoswiss.com)

Web: [www.cryoswiss.com](http://www.cryoswiss.com)

small – easy to handle - mobil

This revolutionary and patented cryosurgical system provides reproducible treatment results for a variety of every day dermatological indications as warts, ceratosis, age pigmentations e.g. by a constant temperature of  $-89\text{ }^{\circ}\text{C}$  (liquid N<sub>2</sub>O) .  
 Within about 10 – 12 sec. The treated skin is frozen in a depth of max. 3 mm. CRYOALFA can be dosed exactly and it is ready for operation anytime and everywhere because of its mobility.



### More advantages

- no destruction of biological functions and structures
- satisfying cosmetic effects
- nearly no scars
- bloodless treatment without anesthesia
- almost painless and therefore also particularly suitable for kids
- good alternative in case of high operation risk or impossibility
- most everyday indications disappearing mostly after only one application

**Super C**



**Super**



**Lux C**



**Lux**



**Perfect C**



**Perfect**



A number of special applicators enlarge the possibilities of applications in the dermatology and other medical spheres.

We are the inventors of the Liquid Freezing of system. Our innovation has proven itself in the practice for more than 15 years around the world.

Our prices are most competitive. Come and see for yourselves! We would like to give advice to you personally.

**Telephone +49 351 795 1866**

**E-Mail: [medica@cryoswiss.com](mailto:medica@cryoswiss.com)**

- Dermatology (Cryodestruction of skin lesions)
- Esthetic medicine
- Gynecology
- Oncology
- Proctology
- Dentistry
- Veterinary medicine
- Cosmetics and medical chiropody
- Minimal invasive therapy



Cryoalfa products are manufactured according to the latest technical standards, comply with the requirements of the Medical Devices Regulations and they have CE certification.

